

H3 Electrocardiograph

Display:

6.2" LCD touch screen (resolution:600 × 480)

3 × 4 and 6 × 2 display format, 12-lead waveforms display simultaneously;

Intelligent brightness sensor for automatic screen brightness adjust



COMeN

Recording and Data printing

- 80mm × 20m Roll paper or Z-fold paper
 - Built-in thermal printer / Optional: external printer
 - Recording format:
 - Built-in printer: 3 × 4、3 × 4+1R
 - External printer: 3 × 4、3 × 4+1R、3 × 4+3R、6 × 2、6 × 2+1R、12 × 1、12 × 1+T
 - 10,000 ECG built-in records with binary image format
 - Connectable with FTP server for patient data remote-control
 - Connectable with PC-based ECG management software
 - Unlimited USB data transfer
- Recording Mode
 - a) Default Print: in automatic mode, ECG waveform is printed.
 - b) One-key manual print: By pressing “print” again, rhythm-lead report is printed manually after the default full information report;
 - c) Extendable print: In case arrhythmia occurs, rhythm-lead report is automatically printed after the default full information report.

- 1/3 channel rhythm lead display
- Support real-time sampling and pre-sampling
- Support synchronous sampling and sequential sampling
- ECG review, 600 seconds waveforms review and print for any 12 leads within 10 seconds

- Automatic standby and timing on-off function
- Li-thium battery
- Scanner connectable
- Smart touch screen with slide and multi-finger operate function
- Smart name and fuzzy ID search function



CE 1639



Shenzhen Comen Medical Instruments Co., Ltd.

P/N: EN-H3-2P-20200302-V1.2

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district, Guangming District, Shenzhen, 518106, Guangdong, China.
Tel:+86-755-2640 8879 400-700-9488 Fax:+86-755-2643 1232 Website: en.comen.com E-mail: info@szcomen.com